

User Guide
TRANSIL
Metabolic Bias Kit
TMP-0120-2096

Version 3, Revision 05

HBO consult GmbH
Tel: +49 341 520 44 0
contact@sovicell.com
www.sovicell.com

Contents

1	Quick Protocol	1
2	Background	2
2.1	Applications of TRANSIL Metabolic Bias Assay	4
3	Basic assay principle	4
4	Kit components.....	7
5	Abbreviations.....	7
6	Equipment	8
7	Assay preparation	9
8	Drug candidate preparation	9
9	Replicates.....	10
10	Assay procedure.....	10
10.1	Compound addition	10
10.2	Incubation.....	10
10.3	Sampling of supernatant	11
11	Sample quantification	11
12	Data analysis	11
12.1	Assay parameters	11
12.2	Compound information	13
12.3	Raw data from sample quantification	13
12.4	Results	15
12.4.1	Membrane affinity	15
12.4.2	Fraction Unbound (fu).....	16
12.4.3	Influence of the $f_u(\text{mic})$ on clearance estimates	17
12.4.4	Detailed measurement results –membrane affinity.....	18
12.5	TRANSIL Quality Index	20

12.5.1	Model fit (intercept)	20
12.5.2	Match of measured versus predicted reference signal (ref).....	21
12.5.3	Correlation coefficient (r^2)	22
12.5.4	Number of outliers or missing data (DP).....	23
12.5.5	Data consistency (C).....	23
12.5.6	Slopes of binding.....	24
12.5.7	Reference treatment.....	24
12.5.8	TRANSIL peak areas exceed reference peak areas.....	24
13	Storage and shelf life	26
14	Trouble shooting.....	26
14.1	Poor recovery	26
14.1.1	Challenges and problem identification	26
14.1.2	Problem-solving approaches.....	26
14.2	Non-linearity of the response.....	27
14.2.1	Challenges and problem identification	27
14.2.2	Problem-solving approaches.....	27
14.3	Low Membrane affinity	28
14.3.1	Challenges and problem identification	28
14.3.2	Problem-solving approaches.....	29
15	References	30

1 Quick Protocol

1. Plate Thawing and preparation

- Thaw plate or individual tube units for 3h at room temperature (alternatively overnight).
- Spin plate quickly for 5 seconds at 750 g.
- Make sure the plate has a working temperature between 20°C and 25°C when starting the experiment.
- Leave caps closed while preparing the test compound.

2. Drug Candidate Preparation

- Prepare 16x stock solutions of each compound in 32% solvent (e.g. DMSO) - yields a final solvent conc. of 2%.
- The final compound concentration in the assay depends on the compounds solubility, analytical method and instrumentation: If permitted by compound solubility use 1 μM final assay concentration. This requires 16 μM stock solutions.
- Since each compound is added in an aliquot of 15 μl to each well of an 8-well tube unit, at least 120 μl stock solution are required for each compound. Allow an additional 80 μl for accurate pipetting.

3. Drug Candidate Addition

- Open wells with supplied decapper.
- Mix the stock solutions carefully.
- Transfer 15 μl of the 16x stock solution to a column of 8 wells of the TRANSIL assay plate proceeding column by column. Change tips after each transfer step to avoid carryover of beads.
- Close tube wells and make sure that the capband is oriented in the same direction as before.

4. Incubation and Supernatant Sampling

- Incubate the plates on a plate shaker at 1000 rpm for 12 minutes.
- Spin the plate in a swing-out centrifuge for 10 minutes at 750 g.
- Transfer 50 – 100 μl supernatant in a standard 96 well plate for analytical quantification. Make sure that no beads are carried along.

5. Analysis

- Quantify supernatants by the method of choice.
- For evaluation of the results, please use the supplied MS Excel spreadsheet and refer to the operating instructions for data analysis.

2 Background

The complex, costly, and often uncertain outcome of the drug discovery and development process requires the simultaneous optimization of several properties. It has now long been recognized that favorable potency and selectivity characteristics are not the sole hallmarks of a successful drug discovery program, nor is the safety profile considered to be the only hurdle to be overcome, although it is of paramount importance.

The ability to prospectively predict the pharmacokinetics of new chemical entities in humans is a powerful means by which scientists involved in the discovery of new drugs can select for further development only those compounds with the potential to be successful therapeutic agents.

The half-life of a drug is a major contributor to the dosing regimen, and it is a function of the clearance and apparent volume of distribution (VD), each of which can be predicted and combined to predict the half-life. Drugs with short half-lives are more likely to be required to be administered more frequently than those with long half-lives. Dosing regimen is also intrinsically linked to other factors such as the pharmacodynamics of the drug and the difference between systemic concentrations associated with side effects vs those minimally required for efficacy. However, these latter attributes are much more difficult to predict from in vitro or animal data and will be different for each therapeutic target. Thus, a great deal of focus has been placed on the prediction of human half-life. While methods using allometric scaling or correlative methods exist for prediction of half-life, greater success is attained if the two major components of half-life, clearance and volume of distribution (note that the TRANSIL Intestinal Absorption Kit is well suited for VD predictions), are predicted separately and combined to generate a half-life prediction (Obach et al 1997).

When estimating intrinsic clearance in metabolic stability incubations with microsomes or hepatocytes, it is often assumed that the total test item concentration is freely available for enzymatic metabolism. However, this is not the case for many compounds. A significant fraction may bind to intracellular membranes – including endoplasmic reticulum, mitochondrial, or plasma membranes – thereby reducing the unbound concentration in the aqueous phase. Since only the unbound drug is available to interact with metabolic enzymes, such binding reduces the effective exposure of the drug to metabolism.

Consequently, it is important to estimate microsomal or hepatocyte binding in these in vitro incubations. For example, if the observed half-life yields a Clint estimate of 5 $\mu\text{L}/\text{min}/10^6$ cells, but only 1% of the drug is unbound ($f_u = 0.01$), the corrected Clint is actually 500

$\mu\text{L}/\text{min}/10^6$ cells — a 100-fold difference. This bias can lead to incorrect ranking of compounds, poor in vitro–in vivo extrapolation, and misjudgment of metabolic stability.

The consequence is a systematic underestimation of intrinsic clearance (CL_{int}) when membrane binding is ignored. The observed metabolic rate is interpreted as if it resulted from the total compound concentration, although in reality, the metabolic turnover was driven only by the unbound fraction. The true intrinsic clearance should be corrected by dividing the observed clearance by the unbound fraction (f_u), using the equation:

$$CL_{int}^* = \frac{CL_{observed}}{f_u}$$

Consequently, it is important to quantify the extent of membrane binding in these in vitro incubations. Traditionally, this has been done by estimating the free fraction of drug using dialysis systems with microsomes (without cofactors). However, dialysis-based methods are time-consuming and often require overnight equilibration. During this prolonged incubation, test items may degrade through hydrolysis or precipitate due to limited solubility — particularly for lipophilic or unstable compounds.

In hepatocyte-based dialysis experiments, an additional challenge arises: metabolic activity may not be fully suppressed. While inhibitors are sometimes added to block enzymatic turnover, the diversity of metabolic enzymes in hepatocytes and the likelihood of residual enzymatic activity make it difficult to ensure complete inhibition. As a result, compound degradation during dialysis may confound the measurement of binding.

Moreover, the free fraction determined by dialysis is specific to the microsomal or hepatocyte protein concentration used in the setup, which often differs from that applied in metabolic stability assays. This mismatch can introduce further error into clearance estimates if not properly accounted for.

Membrane binding reduces not only the concentration of a drug available for metabolism by CYP enzymes, but also the concentration available to inhibit those enzymes. This can result in a substantial underestimation of inhibitor potency, often seen as inflated IC_{50} or K_i values. This effect is particularly pronounced for lipophilic and basic compounds that show high nonspecific binding. As a consequence, the true risk of drug–drug interactions may be overlooked if membrane binding is not considered.

This limitation applies to both microsomal and hepatocyte-based inhibition studies. In mechanism-based inhibition experiments using microsomes, high protein concentrations are commonly used, which increases the extent of membrane binding and reduces the unbound drug concentration. In hepatocyte assays, the situation is further complicated by the diversity

and abundance of intracellular membranes. These additional binding sites reduce the effective intracellular concentration of the test compound, lowering its interaction with metabolic enzymes.

Correcting for the unbound fraction in both microsomal and hepatocyte incubations is therefore essential to accurately evaluate enzyme inhibition and assess the potential for drug–drug interactions.

2.1 Applications of TRANSIL Metabolic Bias Assay

The TRANSIL Metabolic Bias Assay determines the affinity of drug compounds for biological membranes, such as those found in microsomes or hepatocytes. The measured parameter is a distribution coefficient, defined as the ratio of the drug concentration in the lipid phase to that in the buffer phase. This coefficient reflects a fundamental physicochemical property of each compound and enables the calculation of unbound drug fractions across a wide range of experimental conditions.

Because the distribution coefficient is independent of the specific assay microsomal (protein) concentration or hepatocyte density, it can be used to correct for membrane binding in both microsomal and hepatocyte incubations, and to estimate intracellular free concentrations in liver tissue if membrane density is known.

The assay is designed for rapid execution, requiring only 12 minutes of incubation. It is available in a ready-to-use 96-well plate format and is fully compatible with liquid handling systems for automation. This makes it ideally suited for high-throughput applications in drug discovery and ADME profiling.

3 Basic assay principle

The principle of the TRANSIL Metabolic Bias Assay is to assess the affinity of test compounds to biological membranes, including those found in microsomes or hepatocytes (Figure 1). The membrane affinity is determined by incubating a fixed concentration of the drug candidate with varying concentrations of membrane surface area immobilized on the silica beads. A total of 8 wells of a tube unit/plate are used to determine the microsomal membrane affinity for each compound (Figure 2). Six wells contain microsomal membrane silica beads while two serve as references to account for non-specific binding and contain buffer only. Using the spreadsheet and algorithms supplied with the assay, the affinity to the microsomal membranes is calculated from remaining free compound concentration in the supernatant of each well with membrane beads. Any of the available detection systems, such as HPLC-UV,

LC-MS/MS, scintillation counting, etc. can be used for quantification, as long as it can quantify μM concentrations in volumes of 50 μl or less.






Figure 1: TRANSIL Metabolic Bias Assay principle. a) Illustration of the assay workflow: The same amount of drug is added to all wells followed by a mixing step. After 12 minutes incubation, beads are separated from the buffer by centrifugation and the remaining supernatant is sampled for quantification. b) Illustration of a typical metabolic stability assay setup (left panel) compared to a TRANSIL incubation (right panel) containing beads with immobilized microsomal membranes (blue spheres associated with membranes). Only the proportion of compound freely available in the buffer phase can be metabolized by liver microsomal enzymes (metabolism enzymes are illustrated by the red spherical objects located in the blue membrane vesicles). Depending on the affinity of the test compound to microsomal membranes a certain proportion of the compound is distributed within the membrane phase and thus, lowering the free concentration. The $f_u(\text{mic})$ obtained from the TRANSIL experiment can be used to correct the metabolism rates in the stability setup.

a)



b)



c)




Figure 2: The TRANSIL Metabolic Bias Kit uses a column of 8 wells to determine the affinity to microsomal membranes (MIC). a) Photograph of the assay plate and b) the annotated tube units supplied. c) Illustration of the assay plate showing the reference rows A and H (white wells) as well as the increasing microsomal membrane concentration from wells B to G (blue). The dashed line indicates the row orientation of the plate: the same amount of drug is added to all tube wells in one column. The plate can be used for 12 compounds.

4 Kit components

A TRANSIL Metabolic Bias Kit is composed of the following items:

No.	Qty.	Item
1	1	One 96-well plate containing twelve units of 8 tubes, each filled with a suspension of TRANSIL Metabolic Bias Beads in 10 mM phosphate-buffered saline (pH 7.4). The tube units are securely locked into the assay plate to support reliable handling with automated liquid handlers. They can be easily released from the underside of the plate, allowing users to run fewer than twelve test compounds per experiment if desired.
2	1	Decapper-8
3	1	Instruction manual
4	1	CD with spreadsheet calculation


5 Abbreviations

cmp	Compound
conc	Concentration
DMSO	Dimethyl sulfoxide
$f_u(\text{mic})$	Unbound fraction of drug in metabolism experiment
MA	Membrane affinity defined as the concentration of drug in membrane (lipid) over concentration of drug in buffer: $MA = \frac{c_l}{c_b}$. The mass balance equation is used to calculate membrane affinity from experimental data.
PBS	Dulbecco's Phosphate buffered saline used in 1x concentration
TQI	TRANSIL Quality Index
r^2	Correlation coefficient

V_b	Buffer volume
V_l	Lipid volume

6 Equipment

The following equipment is required to run the TRANSIL Metabolic Bias kit:

No.	Instrument	Specification
1	Plate shaker	For high speed mixing (min. 800 rpm), i.e. MixMate (Eppendorf).
		 <p>Alternatively, a vortexer with a plateholder can be used.</p>
2	Centrifuge	Including rotor for SBS standard assay plates

7 Assay preparation

Upon receipt the kit should be stored at -20°C (-4°F).

Before use, thaw the assay at 4°C for a period of 12 hours (overnight) or, at room temperature for a period of 3 hours. Make sure the tubes have reached room temperature (between 20° and 25°C) prior to starting the assay. After thawing, spin plate quickly for 5 seconds at 750 g to collect all liquid at the bottom.

If it is desired to analyze less than 12 compounds at the same time, it is possible to remove columns of 8 tubes, interlocked by the lid-strip. We advise to remove the strips which shall be saved for future experiments and leave the tubes for current use on the rack. Remove tube strips by carefully pushing the individual tubes up from the bottom of the plate rack. Always keep lids closed when removing tubes.

8 Drug candidate preparation

Prepare a 16x stock solution for each drug candidate in DMSO. The final assay DMSO concentration can range from 2% to 6%. A 2% DMSO concentration is recommended (requires 32% DMSO in 16x compound stock) as higher DMSO concentrations may result in slight underestimation of binding.

Please consider the following:

Concentration: The TRANSIL Metabolic Bias Kit can be used in conjunction with different analytical methods and instruments. These include LC-MS/MS, as well as other methods such as scintillation counting. Please note that the lower limit of the compound concentration in the assay is only limited by the detection limit and dynamic range of the analytical system used. However, we advise to choose a compound concentration high enough to assure that the quantification is fully within the linear range of the instrument. Alternatively, it is advised to prepare a detailed calibration curve to account for non-linearities. Please contact the customer service for further advice on the best approach to the particular compound and situation.

The upper limit of the compound concentration in the assay is limited by the compounds solubility as well as the saturation of individual beads or the entire bead suspension with the test compound. Therefore, we recommend using final assay concentrations of 1 µM or less.

Volume: We recommend preparing a volume of at least 200 μl per compound. It is necessary to have at least 120 μl of the stock solution for each compound drug candidate since to each of the 8 tube wells 15 μl of the compound is added.

9 Replicates

The TRANSIL Metabolic Bias assay is designed such that one compound utilizes 8 wells – two references and 6 wells with increasing immobilized biological phase (membrane surface area). Therefore, the assay provides 6-fold determination of the assay parameters. Thus, it is not necessary to run more than one row per compound to obtain replicates for statistical validity.

10 Assay procedure

Follow the following 5 steps for the assay procedure:

10.1 Compound addition

Mix the compound stock solution carefully by vortexing. When the TRANSIL Metabolic Bias kit has reached room temperature and the plate has been centrifuged briefly, remove the capbands with the decapper only immediately before compound addition. Make sure to maintain the original capband direction so that lids will be returned to the original wells to avoid any cross-contamination of beads etc. Add 15 μl of test compound to each well of a tube unit of 8 wells. Use one tube unit per compound (for example wells A1 to H1) so that twelve compounds can be analyzed using one kit. Change tips after each compound transfer step to avoid carryover of beads.

10.2 Incubation

Incubate the plate with repeated aspiration and suspension to ensure proper mixing:

- Small molecules: 20 cycles
- Large molecules: 120 cycles

Alternatively, a plate shaker may be used @ 1000 RPM for 12 minutes at RT. The first time a plate shaker is used for TRANSIL assays it is essential to determine that all the beads are resuspended in solution. To ensure beads are resuspended, visually inspect the plate after 1 Separation of beads and buffer

Spin the plate for 10 minutes at up to 750 g to sediment the beads from the suspension.

10.3 Sampling of supernatant

Take 50 – 100µl samples from the supernatants for analysis. Handling tips:

- Make sure that no beads are carried along when transferring the supernatant to the quantification plate.
- For supernatant sampling we advise not to remove the tubes from the rack. However, it may be convenient to remove and discard closed tube strips after supernatant sampling for easier access to the remaining tubes on the rack. Make sure to close the tubes after sampling and before discarding.
- When manually sampling supernatants we advise to guide the pipette tips along the tube walls.

11 Sample quantification

Use your analytical technique of choice for quantifying the compound concentration in the supernatant obtained in the last assay step.

12 Data analysis

Open the supplied spreadsheet for data analysis and follow the steps below to obtain the results for the TRANSIL Metabolic Bias kit. Only the fields marked in green require user input. Cells marked with gray background contain default values which may need to be adjusted (Figure 3).

12.1 Assay parameters

Open the “main” tab and enter the assay parameters in the column C8 to C12. Enter the lot specific parameters from the certificate of analysis that came with the assay plate. Also, enter the lot number and the assay date.

	A	B	C	D	E	F	G	H	I
1	Transil Metabolic Bias Kit								
2	Version 3 Revision 06								
3			required input						
4			optional input						
5									
6	Assay Plate Parameters		T1:	T2:	T3:	T4:	T5:	T6:	
7			0,029 µl	0,052 µl	0,094 µl	0,170 µl	0,306 µl	0,550 µl	
8									
9		Lot Number:	XXXX						
10									
11		Date:	07.06.24						
12									
13									
14									
15		Compound name	sample volume [µL]	stock conc. [µM]	logMA_{mic}		r²	TQI¹	
16			15,0 µl	80,0 µM					
17	No 1	cmp 1	15,0 µl	80,0 µM	2,43	± 0,05	0,9884	●	9,5
18	No 2	cmp 2	15,0 µl	80,0 µM	2,92	± 0,02	0,9986	●	9,7
19	No 3	cmp 3	15,0 µl	80,0 µM	2,28	± 0,01	0,9969	●	9,7
20	No 4	cmp 4	15,0 µl	80,0 µM	3,74	± 0,05	0,9919	●	7,8
21	No 5	cmp 5	15,0 µl	80,0 µM	4,57	± 0,07	0,9769	●	7,6
22	No 6	cmp 6	15,0 µl	80,0 µM	4,34	± 0,05	0,9927	●	7,8
23	No 7	cmp 7	15,0 µl	80,0 µM	4,17	± 0,05	0,9901	●	7,8
24	No 8	cmp 8	15,0 µl	80,0 µM	3,97	± 0,09	0,9923	●	7,8
25	No 9	cmp 9	15,0 µl	80,0 µM	5,04	± 0,03	0,9737	●	7,7
26	No 10	cmp 10	15,0 µl	80,0 µM	4,53	± 0,04	0,9860	●	7,7
27	No 11	cmp 11	15,0 µl	80,0 µM	3,40	± 0,04	0,9960	●	8,7
28	No 12	cmp 12	15,0 µl	80,0 µM	4,08	± 0,05	0,9904	●	7,8
29	¹ TQI (Transil Quality Index): between 7 and 10: good data quality between 5 and 7: compromised data quality below 5: poor data quality								
30									
31									
32	ID	test article	fu mic (0,5 mg/ml)	fu Hepatocytes					
33				Human	Dog	Rat	Mouse		
34	1	cmp 1	92,8%	99,9%	99,7%	99,7%	99,6%		
35	2	cmp 2	80,6%	99,6%	99,0%	99,2%	98,7%		
36	3	cmp 3	94,8%	99,9%	99,8%	99,8%	99,7%		
37	4	cmp 4	38,5%	97,1%	93,6%	94,9%	91,8%		
38	5	cmp 5	8,4%	83,2%	68,2%	73,1%	62,1%		
39	6	cmp 6	13,7%	89,5%	78,7%	82,4%	73,8%		
40	7	cmp 7	19,0%	92,6%	84,5%	87,3%	80,6%		
41	8	cmp 8	27,1%	95,2%	89,7%	91,7%	86,8%		
42	9	cmp 9	3,1%	63,0%	42,4%	48,3%	36,0%		
43	10	cmp 10	9,2%	84,5%	70,3%	75,0%	64,3%		
44	11	cmp 11	57,9%	98,7%	97,0%	97,6%	96,1%		
45	12	cmp 12	22,4%	94,0%	87,1%	89,5%	83,7%		

Figure 3: Screen shot of the “main” tab of the spreadsheet for analyzing data from the TRANSIL Metabolic Bias Kit. The “main” tab is used to enter lot specific data as well as for reporting final results, the TRANSIL Quality Index (TQI) and predictions of the unbound fraction in microsomal incubations.

12.2 Compound information

Please enter the compound names in the column B17 to B28 of the “main” tab. Enter the sample volume added to each well in field C16. If a different sample volume was used for each drug, enter the sample specific volumes in the column C17 to C28.

Enter the concentration of the samples’ stock solutions in field D16. If a different sample concentration was used for each drug, enter the sample specific concentrations in the column D17 to D28 (remember this is the concentration of the stock solution).

12.3 Raw data from sample quantification

Open the tab “raw data” and enter the peak areas or heights for each well in column G (Figure 4). Note that column A lists the name of the compound used in each well. Caution: Make sure to begin data entry in field G6 for the first well of the plate (A1). When scrolling through the spreadsheet the header line in row 5 remains in place, while the fields for peak area entry move up.

	A	B	C	D	E	F	G	H	I	J	
1	Please enter the peak area or concentration data in column G below										
2											
3	Please leave missing data fields blank										
4											
5	test article		Well	Sample	Area / height	nm / amu	Note				
6	cmp 1		A-1	Ref 1	199880	278.4 / 121.1					
7	cmp 1		B-1	Well 1	202710						
8	cmp 1		C-1	Well 2	193380						
9	cmp 1		D-1	Well 3	184270						
10	cmp 1		E-1	Well 4	166290						
11	cmp 1		F-1	Well 5	155210						
12	cmp 1		G-1	Well 6	127620						
13	cmp 1		H-1	Ref 2	213680						
14	cmp 2		A-2	Ref 1	157880	399.1 / 119.1					
15	cmp 2		B-2	Well 1	150350						
16	cmp 2		C-2	Well 2	143490						
17	cmp 2		D-2	Well 3	131270						
18	cmp 2		E-2	Well 4	108580						
19	cmp 2		F-2	Well 5	83639						
20	cmp 2		G-2	Well 6	58063						
21	cmp 2		H-2	Ref 2	179420						
22	cmp 3		A-3	Ref 1	304310	837.6 / 158.1					
23	cmp 3		B-3	Well 1	318260						
24	cmp 3		C-3	Well 2	314560						
25	cmp 3		D-3	Well 3	303540						
26	cmp 3		E-3	Well 4	289720						
27	cmp 3		F-3	Well 5	267220						
28	cmp 3		G-3	Well 6	227120						
29	cmp 3		H-3	Ref 2	346400						
30	cmp 4		A-4	Ref 1	257950	329.3 / 162.1					
31	cmp 4		B-4	Well 1	101910						
32	cmp 4		C-4	Well 2	76804						
33	cmp 4		D-4	Well 3	56021						
34	cmp 4		E-4	Well 4	37631						
35	cmp 4		F-4	Well 5	28307						
36	cmp 4		G-4	Well 6	17816						
37	cmp 4		H-4	Ref 2	269860						
38	cmp 5		A-5	Ref 1	19699	285.2 / 152.2					
39	cmp 5		B-5	Well 1	1932.1						
40	cmp 5		C-5	Well 2	1332.5						
41	cmp 5		D-5	Well 3	849.27						
42	cmp 5		E-5	Well 4	613.04						
43	cmp 5		F-5	Well 5	455.1						
44	cmp 5		G-5	Well 6	306.7						
45	cmp 5		H-5	Ref 2	30682						
46	cmp 6		A-6	Ref 1	155480	267.2 / 193.3					
47	cmp 6		B-6	Well 1	22414						
48	cmp 6		C-6	Well 2	15924						
49	cmp 6		D-6	Well 3	10316						
50	cmp 6		E-6	Well 4	6649.4						
51	cmp 6		F-6	Well 5	4452.1						
52	cmp 6		G-6	Well 6	2859.6						
53	cmp 6		H-6	Ref 2	156490						
54	cmp 7		A-7	Ref 1	58074	319.3 / 200.2					
55	cmp 7		B-7	Well 1	11758						

Figure 4: Screen shot of the “rawdata” tab of the spreadsheet for analyzing data from the TRANSIL Metabolic Bias Kit. The “rawdata” tab is used to enter peak area or concentration data from the supernatants of the assay plate after incubation and centrifugation.

12.4 Results

The spreadsheet calculates membrane affinities and QC parameters immediately after entering the lot specific information, compound names and concentrations, as well as the raw data from quantification.

12.4.1 Membrane affinity

The membrane affinity is a partitioning coefficient of drug between membrane and buffer. It is defined as the concentration of drug in membrane over the concentration of drug in buffer:

$$MA = \frac{c_1}{c_b} \quad (1)$$

The membrane affinity is calculated from the assay data using the mass balance equation:

$$n_t = c_b \cdot V_b + c_1 \cdot V_1 \quad (2)$$

which is rearranged such that the membrane affinity can be determined from the slope of plotting the ratio of total amount of drug (n_t) over remaining concentration in supernatant (c_b) against the lipid membrane volume present in each well:

$$\frac{n_t}{c_b} = \frac{c_1}{c_b} \cdot V_1 + V_b = MA \cdot V_1 + V_b \quad (3)$$

The results for the membrane affinity are reported in column E17 to E28 along with the TRANSIL Quality Index.

Results with an index greater than 7 are of good quality, results with an index between 5 and 7 are compromised, but may be reasonably accurate, while results with an index below 5 are poor and should be reported with caution.

The default requirement for good references is that both measurements will not deviate more than 30%. This assumption can be changed by setting the margin in cell I11 to a different value. If the references differ more than this threshold of 30% the spreadsheet uses the highest reference value. However, if the highest reference value is lower than the concentration determined in the TRANSIL well with the lowest membrane surface area, then the spreadsheet discards the reference measurements and selects the first TRANSIL measurement as reference and eliminates this TRANSIL measurement from the calculation of

the membrane affinity. When this approach is used, the reported membrane affinity will be higher or equal the true membrane affinity. Please refer to the trouble shooting section if this occurs.

12.4.2 Fraction Unbound (f_u)

The fraction of a test item that remains unbound in an incubation — referred to as f_u — depends on three key parameters: the total membrane volume present in the assay (V_l), the incubation volume (V_{inc}), and the compound's affinity for the membrane (MA), as determined with the TRANSIL Metabolic Bias Kit.

Membrane affinity is defined as the partition coefficient between the membrane phase and the buffer phase (see above). Given this affinity, the unbound fraction can be calculated using the following equation:

$$f_u = \frac{1}{MA \cdot \frac{V_l}{V_{inc} + V_l} + 1}$$

This equation describes how much of the compound remains in the aqueous phase and is thus available to interact with metabolic enzymes. The higher the membrane affinity or the larger the membrane volume relative to the incubation volume, the lower the unbound fraction will be.

The V_l/V_{inc} ratio can be derived from known assay parameters. For microsomal incubations, the membrane volume is typically estimated from the protein concentration using a lipid-to-protein ratio (CLP). For hepatocyte incubations, we have calculated the membrane content from the size and shape and other known parameters of human, dog, rat and mouse hepatocytes. These parameters are preconfigured in the analysis spreadsheet and can be adjusted to reflect the specific conditions of the metabolic stability experiment.

The lipid to protein concentration ratio CLP has been experimentally determined for human liver microsomes:

$$CLP_{\text{human mic}} = \frac{c_l}{c_p} = 0.581 \frac{\mu\text{l}_{\text{lipid}}}{\text{mg}_{\text{protein}}}$$

This ratio is used to calculate the lipid content in the human microsomal incubations given the assay volume V_t and the protein content c_p of the microsomes:

$$V_l = CLP \cdot c_p \cdot V_t$$

In the supplied spreadsheet, the calculation of f_u is automated based on assay-specific inputs. For microsomal incubations, the user should enter the microsomal protein concentration in cell F32, and check whether the default lipid-to-protein ratio (CLP) of 0.5 mg/mL in cell L32 is applicable to their batch. For hepatocyte incubations, the user is required to enter the cell density (in cells/mL) in cell N32, which will be used together with a default hepatocyte CLP to estimate the intracellular membrane content.

Microsomal incubations can vary significantly in lipid content depending on the species and even from lot to lot within the same species. This variability affects the lipid-to-protein ratio (CLP), which in turn influences the calculated membrane volume and the resulting unbound fraction. The CLP can be determined by straightforward lipid quantification procedures. We offer a standard protocol for lipid analysis upon request. Alternatively, lipid content analysis of microsomal or cellular suspensions is available as a service.

12.4.3 Influence of the $f_u(\text{mic})$ on clearance estimates

Since only the unbound fraction of drug is available for metabolic decay the clearance estimated in the microsomal incubation Cl_{inc} is lower than the theoretical intrinsic clearance Cl_{int} that would be observed if the drug would not bind to microsomal membranes. The equation

$$Cl_{int} = \frac{Cl_{inc}}{f_u}$$

describes the relationship between the clearance rate measured in the incubation and the intrinsic clearance. However, one should note that the in vivo clearance will also be influenced by reduction of the free fraction of drug in the liver tissue due to binding to plasma membranes. The most significant membranes are the following:

1. Endoplasmic Reticulum (ER) (60–70%)

The most extensive membrane system in hepatocytes, especially the smooth ER, which houses cytochrome P450 enzymes. It is the principal contributor to nonspecific drug binding in metabolic assays.

2. Mitochondrial Membranes (15–20%)

Includes both outer and inner mitochondrial membranes. These contribute substantially to total intracellular surface area and are relevant binding sites for lipophilic compounds.

3. Plasma Membrane (5–10%)

While structurally external to the cell, the plasma membrane significantly contributes to overall drug partitioning. It plays a role in both initial uptake and intracellular distribution.

4. Golgi Apparatus (3–5%)

A stacked membrane system involved in protein processing and transport. It adds a moderate amount of surface area to the intracellular membrane pool.

5. Lysosomal Membranes (2–3%)

6. Endosomal Membranes (1–2%)

7. Peroxisomal Membranes (<1%)

8. Nuclear Envelope (<1%)

Note that percentages are approximate and may overlap depending on structural definitions and measurement methods. The total may slightly exceed 100% due to classification overlap.

12.4.4 Detailed measurement results –membrane affinity

Detailed measurement results can be found for each drug can be found on the spreadsheet's detail tabs with the indices from 1 to 12 for each respective drug. Figure 5 illustrates the information reported on each individual drug tested.




Figure 5: Screen shot of the details “1” tab of the spreadsheet for analyzing data from the TRANSIL Metabolic Bias Assay kit. The “1” tab shows calculated concentrations in each well and all calculations performed to derive the affinity to membranes as well as three plots indicating the experiment performance.

12.5 TRANSIL Quality Index

The TRANSIL Quality Index (TQI) is based on five independent measures derived from the data analysis. For each individual measure a partial quality score on a scale between 0 and 10 is attributed to the estimate. 0 represents lowest quality, while 10 represents highest quality. The final quality index is a weighted average of the partial quality scores.

12.5.1 Model fit (intercept)

The membrane affinity is calculated by fitting the experimental data to the rearranged mass balance equation:

$$\frac{n_t}{c_b} = MA \cdot V_l + V_b \quad (3)$$

Fitting optimal data to equation (3) will yield a slope that exactly represents the true membrane affinity, MA, and the buffer volume used in the experiment. In fact, a biased estimation of the slope will typically result in a biased estimation of the intercept as well. Since the intercept equals the buffer volume used in the experiment, the estimated intercept is used as a quality control parameter. If the estimated buffer volume is within an interval $\pm 10\%$ around the true value a partial quality score of 10 is attributed. If the estimated buffer volume is within an interval $\pm 50\%$ around the true value a partial quality score of 5 is attributed. The partial quality score for the model fit has a weight of 3 in the total quality index.




Figure 6: Illustration of fitting experimental data to equation (3) to determine the membrane affinity. A buffer volume of 240 μL has been used in the experiment, the intercept was estimated to 346 μL , hence a quality score of 5 was attributed to the model fit.

12.5.2 Match of measured versus predicted reference signal (ref)

When determining the membrane affinity via the six different lipid volumes using TRANSIL beads along with 2 reference estimates without TRANSIL beads, the expected peak area resulting from quantification of the references can be calculated from the peak areas from the TRANSIL wells by linear regression, since lipid binding can be assumed to be a non-cooperative process (Figure 7). This score has a weight of 3 in the TQI.




Figure 7: Illustration of estimating the peak area in the reference vials by plotting the inverse of the peak areas of the compound concentration of supernatants in TRANSIL vials against the lipid content. The inverse of the intercept represents the expected peak area of the references.

A deviation of the measured from the expected reference peak area can be due to a non-linear calibration curve or unspecific binding, which will be more pronounced in the references without the lipid phase of the TRANSIL beads than in the TRANSIL wells. Table 1 lists the partial quality scores for deviations of the reference peak areas from expected reference peak areas.

Table 1: Partial quality scores for deviations of the reference peak areas from expected reference peak areas.

Deviation	Score
10.0%	10
20.0%	9
50.0%	8
100.0%	7
200.0%	5
500.0%	3
>500.0%	0

12.5.3 Correlation coefficient (r^2)

The correlation coefficient from fitting the experimental data to equation (3) also contributes as a partial quality score (Table 2). This score has a weight of 3 in the TQI.

Table 2: Partial quality scores for the least square model fit of the experimental data to equation (3).

r^2	Score
0.9999	10
0.999	9
0.99	8
0.9	7
0.8	6
0.7	5
0.6	4
0.5	3
<0.5	0

12.5.4 Number of outliers or missing data (DP)

The number of data points used to calculate the membrane affinity is also used as partial quality score (Table 3). This score has a weight of 2 in the TQI.

Table 3: Partial quality scores for the number of data points used in the model fit of the experimental data to equation (3).

Data points	Score
5	10
4	9
3	6
2	1
1	0

12.5.5 Data consistency (C)

With increasing lipid volume, i.e. increasing lipid membrane surface, the binding of the test items to the membrane should increase proportionally. At least the binding should increase with increasing lipid volume. If the measured peak area suggests decreased binding compared to binding in the TRANSIL well with the next lower lipid volume, then this data point is considered to be inconsistent with the fundamental assumption about lipid binding. If this happens for more than one consecutive TRANSIL well, the data point will be excluded from the calculation. Irrespective of inclusion or exclusion, a partial quality score will be attributed to the data set based on consistency according to Table 4. This score has a weight of only 1 in the TQI as it may affect also the number of data points.

Table 4: Partial quality scores for the number of consistent data points used in the model fit of the experimental data to equation (3).

No. of consistent data points	Score
5	10
4	5
3	2
2	0

12.5.6 *Slopes of binding*

Data fitted to equation (3), plotted in Figure 6 as well as the percentage binding shall increase with increasing membrane surface area (Figure 7). Hence, the slopes of these graphs must all be positive. Most critical of all is the relationship of equation (3), if it has a positive slope it receives a vote of 10 points, otherwise zero. If the relationship plotted in Figure 7 has a positive slope, a vote of 5 points is granted. If the binding curve has a positive slope, a vote of 5 points is granted, otherwise zero. If the total count of votes is 20, a partial score of 10 will be attributed, if the total vote is 10 a score of 5 is attributed, and if the total vote is 0, a partial score of 0 is attributed to the data set.

12.5.7 *Reference treatment*

For each compound two references are measured in the assay kit. If the references vary by no more than 30% and have a higher peak area than the measurements in the TRANSIL wells, the average references is computed and a partial score 10 is attributed. If the reference peak areas are higher than those of the TRANSIL wells, but differences between the two measurements exceed 30%, the maximum of the measurements is chosen. However, if the reference peak areas do not exceed the peak areas from the TRANSIL wells the reference measurements are discarded and the first TRANSIL measurement is taken as reference. In this case a partial score of 6 is attributed.

12.5.8 *TRANSIL peak areas exceed reference peak areas*

The reference peak areas should always exceed the peak areas from the TRANSIL well. If not, the stability or solubility of the compound is compromised. Therefore, the fewer TRANSIL measurements meet this criterion, the lower the partial score attributed to the data set (Table 5).

Table 5: Partial quality scores for the number of data TRANSIL peak areas being higher than reference peak areas.

No. of TRANSIL peak areas higher than reference peak areas	Score
5	10
4	7
3	4
2	2
1	1
0	0

13 Storage and shelf life

The assay kits are shipped in a frozen state and should be stored at -20 °C. TRANSIL materials are stable for several months when stored as recommended. Once thawed and at room temperature, the kit should be used within 24 h.

14 Trouble shooting

14.1 Poor recovery

14.1.1 Challenges and problem identification

Poor data quality such as low TQI's, poor regression fits, or strong variation in duplicate measurements of references may indicate reduced recovery due to poor solubility or stickiness of the test compound. This can result in lower compound concentrations in the reference wells than in the TRANSIL wells. The spreadsheet detects if reference measurements are lower than the signal in the first TRANSIL well. In this case, the spreadsheet replaces the reference value with the measurement from the first TRANSIL well. Consequently, the first TRANSIL well is discarded from the data analysis. Treatment of the references is reported on each compound page in cell D5. The letter "A" (=average) refers to normal treatment as before, "M" (=maximum) is chosen when the difference between references exceeds the value specified in cell I11 of the summary page, and "R" denotes the replacement with the signal in the first TRANSIL well.

For evaluation of recovery issues, include a separate control vial with pure organic solvent (e.g. DMSO) and the test compound in the same concentration as the final assay concentration. Comparison of the peak areas or counts from this organic solvent control and the peak areas from the according calibration signal or the assay references yields a good indication of compound losses through incomplete solubility in the aqueous buffer system or through unspecific binding. Please note that comparing the absolute peak area should be done with caution because of matrix effects.

14.1.2 Problem-solving approaches

- i. Sovicell support team can assist you in checking the plausibility of the data if solubility/non-specific binding problems are observed. In any case, for optimization of the assay parameters it will be helpful to know the solubility of the test compounds in pure buffer solutions.

- ii. DMSO content can be increased. The assay tolerates up to 10% DMSO.
- iii. Test compound concentration can be reduced, however, it has to be considered that running the assay with lower compound concentrations increases the likelihood of measurements outside the linear range of the instruments (c.f. section 15.2).

Before repeating the whole assay you may check the success of recommendations given in ii. to iii. by setting up an individual small control experiment. It is recommended to use the same assay buffer to ensure comparability. Please contact Sovicell support to receive tubes with assay buffer with an appropriate volume.

14.2 Non-linearity of the response

14.2.1 Challenges and problem identification

Frequently, it is observed that mass spectrometers exhibit a non-linear response even in concentration ranges up to 100x above the detection limit. Likewise, impurities of radiolabelled compounds can lead to similar effects when the impurity exhibits different binding properties from the parent compound.

The warning message poor intercept fit or a non-linear shape of the regression (visualized by the “MA Estimation plot” in the individual data analysis tabs of the spreadsheet; see Figure 8) may indicate non-linear response issues.

14.2.2 Problem-solving approaches

- i. Increasing the test compound concentration will increase supernatant concentrations and help to eliminate the non-linear instrument response at low concentrations.
- ii. Non-linear response issue is primarily observed with high affinity compounds. A kit with lower lipid content will increase supernatant concentrations and help to eliminate the non-linear instrument response. To further improve the measurement accuracy of compounds with high membrane affinities we offer the TRANSIL Metabolic Bias Kit for high affinity compounds (TMP-0120-2296).
- iii. If test compound concentration is limited by poor compound solubility, a detailed calibration curve covering the non-linear response can be recorded and used to calculate test compound concentrations. The concentrations calculated from the non-linear calibration curve can then be entered in the spreadsheet’s raw data tab instead of peak areas. Feel free to contact our technical support for guidance, in particular,

because we advise to use the same buffer system for the calibration curve as for the assay.




Figure 8: Illustration of the non-linear response issue which can be recognized by a poor intercept fit (green circles) which usually comes along with a curved plot of $n(\text{total})/c(\text{buffer})$ versus lipid volume plot for the MA estimation. Both the poor intercept fit and the deviation from linearity in this plot are a good indication of the non-linear instrument response to decreasing compound concentrations.

14.3 Low Membrane affinity

If compound binding to TRANSIL is not increasing with increasing TRANSIL concentration, then the compounds exhibit very low affinity to the TRANSIL lipid membrane. This means their membrane affinity is very low. The spreadsheet will automatically use an appropriate alternative approach for the calculation of the membrane affinity if such problems occur.

14.3.1 Challenges and problem identification

Compounds with very low membrane affinity ($\log\text{MA} < 2$) are not accurately measured. Low affinity compounds yield supernatant concentrations in the assay that deviate only marginally from the reference signals (Figure 9).




Figure 9: Illustration of a bar plot of a compound exhibiting low membrane affinity (c.f. individual data analysis tabs of the spreadsheet). The blue bars show the detected signals in the supernatants of TRANSIL wells I to VI. As the compound distributes only weakly into the membranes, supernatant concentrations differ only marginally from the reference signals (red line).

14.3.2 Problem-solving approaches

In case the membrane affinity turns out to be below 300 ($\log MA < 2.5$) we recommend using the TRANSIL Metabolic Bias Kit - low affinity compounds (Product No. TMP-0120-2196). This kit contains an adjusted TRANSIL content (higher than the standard kit) and hence, provides more accurate results.

Technical Support

Phone: +49 341 52044-0

Email: contact@sovicell.com

15 References

- Austin, R. P., Barton, P., Cockcroft, S. L., Wenlock, M. C., and Riley, R. J. (2002). The influence of nonspecific microsomal binding on apparent intrinsic clearance, and its prediction from physicochemical properties. *Drug Metab Dispos* **30**, 1497-1503.
- Hallifax, D., and Houston, J. B. (2006). Binding of drugs to hepatic microsomes: comment and assessment of current prediction methodology with recommendation for improvement. *Drug Metab Dispos* **34**, 724-726; author reply 727.
- McLure, J. A., Miners, J. O., and Birkett, D. J. (2000). Nonspecific binding of drugs to human liver microsomes. *Br J Clin Pharmacol* **49**, 453-461.
- Obach, R. S. (1999). Prediction of human clearance of twenty-nine drugs from hepatic microsomal intrinsic clearance data: An examination of in vitro half-life approach and nonspecific binding to microsomes. *Drug Metab Dispos* **27**, 1350-1359.